

## 9. MAINTENANCE

We recommend having it serviced every 3 to 5 years depending on how often the scales are used.

### ATTENTION!

#### Incorrect measurements due to improper servicing

- Make sure that maintenance and repair are only carried out by an authorised service partner.
- You can find your local service partner at [www.seca.com](http://www.seca.com) or send an e-mail to [service@seca.com](mailto:service@seca.com).

## 10. TECHNICAL DATA

Technical data	
Dimensions <ul style="list-style-type: none"> <li>• Depth</li> <li>• Width</li> <li>• Height</li> </ul>	474 mm 428 mm 2318 mm
Weight	5 kg
Temperature range	+10° C to +40°C
Height of figures	12 mm
Power supply	Battery
Power consumption of head slide <ul style="list-style-type: none"> <li>- with deactivated wireless module and without background lighting</li> <li>- with activated wireless module and permanent background lighting (brightness: 100%)</li> </ul>	20 mA  80 mA
Maximum running time of head slide <ul style="list-style-type: none"> <li>- with deactivated wireless module and without background lighting</li> <li>- with activated wireless module and permanent background lighting (brightness: 100%)</li> </ul>	approx. 3800 Minuten  approx. 2200 Minuten
Measuring data, length measurement <ul style="list-style-type: none"> <li>• Measuring range</li> <li>• Graduation</li> <li>• Accuracy</li> </ul>	30 - 220 cm 1 mm ± 2 mm
Medical device according to Directives 93/42/EEC and 2007/47/EC	Class I with measuring function

Technical data	
Wireless transmission <ul style="list-style-type: none"> <li>• Frequency band</li> <li>• Transmission power</li> <li>• Fulfilled standards</li> </ul>	2.433 MHz - -2.480 MHz <10 mW EN 300328

# 11.ACCESSORIES

Accessories	Article number
Wireless network <b>seca 360° wireless</b> : <ul style="list-style-type: none"> <li>• wireless printer               <ul style="list-style-type: none"> <li>- <b>seca 360° Wireless Printer 465</b></li> <li>- <b>seca 360° Wireless Printer Advanced 466</b></li> </ul> </li> <li>• PC software               <ul style="list-style-type: none"> <li>- <b>seca analytics 105</b></li> </ul> </li> <li>• USB wireless module               <b>seca 360° Wireless USB adapter 456</b> </li> </ul>	country-specific variants country-specific variants Application-specific licensed models 456-00-00-009

# 12.DISPOSAL

## 12.1 Disposal of device



Do not discard with household waste. The device must be disposed of properly as electronic waste. Comply with the national provisions applicable in your country. For further information contact our service department at:

**service@seca.com**

## 12.2 Batteries

Spent (rechargeable) batteries should not be discarded with household waste, regardless of whether they contain harmful substances or not. As a consumer you are obliged by law to dispose of (rechargeable) batteries via the collection points set up by the municipal authorities or the retail sector. Only discard (rechargeable) batteries when fully discharged.